



Fixation Products IV Securement

A complete line of quality fixation products to meet all of your IV securement needs. One company with one goal - to satisfy the customer in the most cost-effective method possible.



STAT 2[®]
CONTROLLER



STAT 2[®] PUMPETTE
CONTROLLER



STAT 2[®] EXTENSION SET
CONTROLLER



To learn more about these and other innovative products, call **800-448-6506** or visit **conmed.com**.

VENI-GARD® IV STABILIZATION DRESSING

FIXATION PRODUCTS



705-4431
Veni-Gard TM



730-4432
Veni-Gard TM Jr.



710-2730
Veni-Gard SP



740-1000
Neoderm



750-4181
Veni-Gard ML



740-1440
Veni-Gard AP



745-1441
Veni-Gard TPN



745-1000
Neoderm

The ultimate in fixation products, providing stabilization, anchoring, and securement in one product.

- Aggressive adhesive foam border system holds objects securely in place and reduces the need for additional securement devices.
- Clear membrane window allows visualization of the site while fixation is maintained.
- Window has a high Moisture Vapor Transmission Rate (MVTR) for breathability.
- Water resistant allowing patients to maintain ADLs.
- Design decreases chance of runaways or catheter pistoning.
- One hand application provides for ease of use.
- Supplied with stretchable sterile tape strips for added stabilization of infusion lines.
- Available in a variety of shapes and sizes for a multitude of fixation applications.
- Sterile.

ClearSite® TM TRANSPARENT MEMBRANE DRESSING

FIXATION DEVICE



A completely transparent fixation device available in a variety of sizes for varying applications.

- Stabilization carrier prevents dressing from curling and self-adhesion upon patient application.
- Unique delivery system allows for true one-handed application.
- High-MVTR ensures breathability and patient comfort.
- Eight sizes, ranging from 2" x 2" to 8" x 10".
- Sterile.

ClearSite® TM	3M Tegaderm™ ¹	S&N OpSite™ ^{1 3}
741-1022	1622W	4654
741-1023	1624W	4007
741-1034	1635	—
741-1044	1630	—
741-1040	1627	4542
741-1045	1626W	4008
741-1068	1628	4631
741-1080	1629	4987

1. ConMed does not manufacture Tegaderm™ or OpSite™ products.

2. Tegaderm™ is a registered trademark of 3M.

3. OpSite™ is a registered trademark of Smith & Nephew.

FIXATION PRODUCTS APPLICATION CHART

	IV Catheter Securement	Peripheral Line Securement	Jugular Line Securement	Central Line Securement	PICC Line Securement	Femoral Line Securement	Epidural Line Securement	Pediatric Line Securement	Fragile Skin Securement	Post-Op Wound Closure	Aggressive Securement
Veni-Gard®											
705-4431		■	■								
730-4432								■			
710-2730	■			■							
750-4181			■	■							■
740-1440					■						
745-1441			■			■	■				
740-1000											■
745-1000					■		■				
ClearSite® TM											
741-1022	■							■			
741-1023	■								■		
741-1034	■	■				■			■		
741-1044		■	■	■	■	■			■		
741-1045				■	■					■	
741-1040										■	
741-1068							■			■	
741-1080										■	

Fixation Products Ordering Information

Veni-Gard®	Description	Metrics	Pouch	Box	Case
701-0510 Strip-N-Tape®	Non-sterile foam tape strip	—	5	1000	5000
705-4431 Veni-Gard® TM	Traditional omega shape	—	1	100	500
730-4432 Veni-Gard® TM Jr.	Pediatric omega shape	—	1	100	500
710-2730 Veni-Gard® SP	Omega shape with IV notch	—	1	100	500
750-4181 Veni-Gard® ML	Elongated omega shape	—	1	25	150
740-1440 Veni-Gard® AP	2" x 3" Narrow bordered	50.8mm x 76.2mm	1	100	500
745-1441 Veni-Gard® TPN	4" x 5" Narrow bordered	102mm x 127mm	1	25	150
740-1000 Neoderm®	2" x 3" Wide bordered	50.8mm x 76.2mm	1	100	500
745-1000 Neoderm® Plus	4" x 5" Wide bordered	102mm x 127mm	1	100	500
ClearSite®	Description	Metrics	Pouch	Box	Case
741-1022 ClearSite® TM	2" x 2" Transparent membrane	50.8mm x 50.8mm	1	100	400
741-1023 ClearSite® TM	2" x 3" Transparent membrane	50.8mm x 76.2mm	1	100	400
741-1034 ClearSite® TM	3" x 4" Transparent membrane	76.2mm x 102mm	1	100	400
741-1044 ClearSite® TM	4" x 4" Transparent membrane	102mm x 102mm	1	100	400
741-1045 ClearSite® TM	4" x 5" Transparent membrane	102mm x 127mm	1	50	200
741-1040 ClearSite® TM	4" x 10" Transparent membrane	102mm x 254mm	1	50	200
741-1068 ClearSite® TM	6" x 8" Transparent membrane	152mm x 203mm	1	25	100
741-1080 ClearSite® TM	8" x 10" Transparent membrane	203mm x 254mm	1	50	200