



Reflex® ELC

Disposable Laparoscopic Clip Applier



DETACHATIP® II



UNIVERSAL PLUS®



ALTRUS®



Learn more about this and other innovative products.

Call 800-448-6506,
visit conmed.com or
email info@mail.conmed.com

Reflex® ELC Laparoscopic Clip Applier


Maximum Visibility

By comparison, the Reflex® ELC jaw design has a 20% greater angle allowing for maximum visibility for precise clip placement.



Maximum Control

The Reflex® ELC has an automatic, direct drive clip feeding system with an audible feedback mechanism indicating complete clip closure. Our superior design provides maximum control, feel and ease-of-use.



Tissue Lock™ Clip

The Tissue Lock™ Clip is designed with opposing male/female interlocking features. This design provides superior on-vessel clip retention and improved hemostasis.

Ordering Information

Reference No.	Product Description	Packaging
530	Reflex® ELC Clip Applier with 20 ML Clips	1pkg, 6/cs