



EndoWeave™



Fenestrated  
(Duckbill)



Allis



Maryland



Tapered  
(Dolphin)



5mm Babcock



Meeker  
(Right Angle)

# LAPAROSCOPIC GRASPERS & DISSECTORS

The High Quality, Cost Effective Alternative



# Laparoscopic Graspers & Dissectors

- 4mm UNIVERSAL Cord eliminates gapping at connection, and the need for other monopolar cords.
- Extended insulation covers the pivot points to reduce the chance of alternative site burns.
- Comfortable, ergonomic handle minimizes fatigue.
- Industry standard 4mm familiar top-post electro-surgical connection
- Easy release trigger eliminates over pressure on release of ratchet
- Finger loop provides extra leverage and comfort
- Available in a wide array of precision tip styles.



Jaw Style	Products	Grasper Ref. # (with ratchet)		Dissector Ref. # (without ratchet)	
		33cm Length	43cm Length	33cm Length	43cm Length
	Tapered (Dolphin)	60-6045-001	60-6045-431	60-6040-001	60-6040-431
	Maryland	60-6045-002	60-6045-432	60-6040-002	60-6040-432
	Fenestrated (Duckbill)	60-6045-003	60-6045-433	60-6040-003	60-6040-433
	Meeker (Right Angle)	60-6045-004	60-6045-434	60-6040-004	60-6040-434
	EndoWeave™	60-6045-005	60-6045-435	60-6040-005	60-6040-435
	Allis	60-6045-006	60-6045-436	60-6040-006	60-6040-436
	5mm Babcock	60-6047-001	60-6047-431	n/a	n/a

- All instruments are packaged sterile, each with a Universal Cord in a peel top tray
- All instruments are Single-Use only
- 1 item/package – 5 packages/case

## CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com

