

**NEW**  
**GoldLine® Pencils**



# GoldLine® Electrosurgical Pencils

# GoldLine® Electrosurgical Pencils

Redesigned with the surgeon in mind.

- Slim, tapered ergonomic design
- Slip-proof handling
- Switch and buttons designed for smooth activation
- Soft blister packaging is puncture and tear resistant; creates less waste
- Available with UltraClean® coated blades preloaded

| Catalog #                   | Product Description                  | Cable Length | Packaging |
|-----------------------------|--------------------------------------|--------------|-----------|
| <b>Hand Switch Pencils</b>  |                                      |              |           |
| Push Button                 |                                      |              |           |
| 131309A                     | UltraClean® Blade, Holster           | 10' Cable    | 40/case   |
| 131326A                     | UltraClean® Blade, Holster           | 15' Cable    | 40/case   |
| 130309A                     | Standard Blade, Holster              | 10' Cable    | 40/case   |
| 130326A                     | Standard Blade, Holster              | 15' Cable    | 40/case   |
| 130303A                     | Standard Blade, Holster, Tip Cleaner | 10' Cable    | 40/case   |
| 130302A                     | Needle Electrode, Holster            | 10' Cable    | 40/case   |
| 138114A                     | Extended Insulated Blade, Holster    | 10' Cable    | 40/case   |
| <b>Rocker Switch</b>        |                                      |              |           |
| 131307A                     | UltraClean® Blade, Holster           | 10' Cable    | 40/case   |
| 131327A                     | UltraClean® Blade, Holster           | 15' Cable    | 40/case   |
| 130307A                     | Standard Blade, Holster              | 10' Cable    | 40/case   |
| 130327A                     | Standard Blade, Holster              | 15' Cable    | 40/case   |
| 130304A                     | Standard Blade, Holster, Tip Cleaner | 10' Cable    | 40/case   |
| 130301A                     | Needle Electrode, Holster            | 10' Cable    | 40/case   |
| 138115A                     | Extended Insulated Blade, Holster    | 10' Cable    | 40/case   |
| <b>Foot Control Pencils</b> |                                      |              |           |
| 138650A                     | UltraClean® Blade, Holster           | 10' Cable    | 40/case   |
| 138652A                     | Needle Electrode, Holster            | 10' Cable    | 40/case   |
| 138659A                     | Blade Electrode, Holster             | 10' Cable    | 40/case   |
| 138655-15A                  | Blade Electrode, No Holster          | 15' Cable    | 40/case   |
| <b>Tip Cleaner</b>          |                                      |              |           |
| 138029                      | Tip Cleaner, Single Use, Sterile     | n/a          | 40/case   |

Please contact your local ConMed Territory Manager for information on our complete product line and accredited training programs.

## CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com

