



Suction Instruments & Tubing



Suction Connecting Tubing

Suction Connecting Tubing serves as the conduit from which suction instruments evacuate fluid and debris into the vacuum canister and source. CONMED's Suction Connecting Tubing is available in Premium Connected, and Economy (100', 30.5m rolls) styles. In addition, a wide range of sizes and configurations allow for efficient and cost-effective setup.

Premium Connected

- Versatile connectors provide a secure fit
- Strong to provide high collapse resistance
- Soft, pliable and low memory for ease of use
- Available with male connectors
- Variety of pre-cut lengths from 1 1/2" (.46m) to 25' (7.6m)
- Three I.D. sizes — 3/16" (4.8mm), 1/4" (6.4mm) & 9/32" (7.1mm)
- Sterile (unless otherwise indicated)

	Ref. No.
3/16" (4.8mm) I.D. x 1 1/2' (.46m) Long, 50/case.....	0036470
3/16" (4.8mm) I.D. x 1 1/2' (.46m) Long, clean, non-sterile, 50/case	0048200
3/16" (4.8mm) I.D. x 3' (.9m) Long, 50/case	0037850
3/16" (4.8mm) I.D. x 6' (1.8m) Long with Male connector, 50/case	0034280
3/16" (4.8mm) I.D. x 6' (1.8m) Long, 50/case	0036480
3/16" (4.8mm) I.D. x 6' (1.8m) Long with Male & Female connectors, 50/case...	0036580
3/16" (4.8mm) I.D. x 10' (3m) Long with Male connector, 50/case	0034310
3/16" (4.8mm) I.D. x 10' (3m) Long, 50/case	0036770
3/16" (4.8mm) I.D. x 12' (3.7m) Long with Male connector, 20/case	0034340
3/16" (4.8mm) I.D. x 12' (3.7m) Long, 20/case	0036560
3/16" (4.8mm) I.D. x 20' (6.1m) Long, 20/case	0036570
1/4" (6.4mm) I.D. x 6' (1.8m) Long with Male connector, 50/case	0034290
1/4" (6.4mm) I.D. x 6' (1.8m) Long, 50/case	0036280
1/4" (6.4mm) I.D. x 6' (1.8m) Long with Male & Female connectors, 50/case...	0037860
1/4" (6.4mm) I.D. x 10' (3m) Long with Male connector, 50/case	0034320
1/4" (6.4mm) I.D. x 10' (3m) Long, 50/case	0036290
1/4" (6.4mm) I.D. x 12' (3.7m) Long with Male connector, 20/case	0034350
1/4" (6.4mm) I.D. x 12' (3.7m) Long, 20/case	0036550
1/4" (6.4mm) I.D. x 20' (6.1m) Long, 20/case	0036540
1/4" (6.4mm) I.D. x 25' (7.6m) Long, 20/case	0036530
9/32" (7.1mm) I.D. x 6' (1.8m) Long with Male connector, 50/case	0034300
9/32" (7.1mm) I.D. x 6' (1.8m) Long, 50/case	0037040
9/32" (7.1mm) I.D. x 9' (2.7m) Long, 20/case	0037050
9/32" (7.1mm) I.D. x 12' (3.7m) Long with Male connector, 20/case	0034360
9/32" (7.1mm) I.D. x 12' (3.7m) Long, 20/case	0037060
9/32" (7.1mm) I.D. x 20' (6.1m) Long, 20/case	0037070



ConMed also offers a broad variety of these products in bulk, non-sterile formats.

The Difference is Clear®