

TITAN 11

All-in-One Digital Radiography

TITAN 11 is a compact digital radiography system which enables various projections possible in any confined operation space. Its configuration will assure smooth transition from analog into digital radiography system while relieving your burden on limited space and budget issue.

Features

- Retrofit System
- Full Automatic operation system
- Remote control
- FPD/CCD detector setting available
- Simple installation



Application



Equipment component



Specification

40/50kW, High Frequency Generator	9M Pixels	C.s.I Scintillator
Maximum 125 kVp	Amorphous Silicon	Manual Collimator
Maximum 500mA	17" x 17" Image size	Synchronization (Up & Down)

Option

- FPD/CCD detector available
- AEC(Auto Exposure Control)
- Carbon table
- DAP
- Auto Collimator

GEMSS MEDICAL SYSTEMS CO., LTD.

61, Dunchon-daero 541 beon-gil, Jungwon-gu,
Seongnam-si, Gyeonggi-do, Korea
Phone: +82-31-764-7321
Fax: +82-31-764-7324
www.gemss-medical.com

The Power of Passion

GEMSS