



# MO1

Pulse Oximeter



# MO1

## Pulse Oximeter



- Streamline design, easthetic & generous.
- Light weight, easy operation.
- 1.8" high definition color screen, Real-time showing of SpO2 waveform, SpO2 & PR value.
- Unique digital SpO2 technology, fast, accurate, stable and strong anti-interference performance.
- At least 48 hours parameter & trend data storage, maximum up to 999 measuring objects.
- It could be connected with PC to download and analyze data.
- Suitable for adult, pediatric and neonate.
- Real-time clock, easy time setup, it could be set manually or via PC software.
- Quick charge via unique seat-charger, real-time battery capacity indicator.
- 10-12 hours battery back-up after complete charge.

## Technical Specifications

### SpO2

Patient type	Adult, pediatric, neonate
Measuring Range	0 ~ 100 %
Resolution	1 %
Accuracy	±2dbm
Parameter	SpO2, PR

### Pulse Rate

Measuring Range	30 ~ 250 bpm
Resolution	1 bpm
Accuracy	±2bpm

### Work Environment

Temperature range operating	0°C ~ 50°C
Transport/Storage	-40° ~ 55°C

### Relative Humidity

Operating	0% ~ 93%
Transport/Storage	0% ~ 95%

### Power Requirements

Operating	100~240VAC, 50/60Hz
Power consumption	<6VA

# Mediblu

3016 NW 82nd Ave  
Doral, FL 33122

Phone: +1 305 406 3931

Fax: +1 305 406 2144

[www.medimarcorp.com](http://www.medimarcorp.com)